

EZRELIEF™

PAIN RELIEF

JUST GOT **EZ'r**

Compression material allows you to stay mobile while recovering

Versatile: Knee, calf, elbow, wrist, thigh, and ankle

Reusable bag maintains product quality

*Product maintains cold skin temperature up to 2X longer than competition

EZRelief™ GEL TECHNOLOGY

UP TO **2x**
LONGER LASTING RELIEF*

Flexible gel that holds temperature longer, helping the body to recover

The diagram shows a cross-section of a gel layer between two fabric layers. A red arrow labeled 'WARM' points down into the gel, and a blue arrow labeled 'COOL' points up into the gel. The gel is depicted as a light blue, textured block.

Hot/cold therapy helps reduce pain and swelling

Additional fabric allows for easy on/off

The Mueller EZRelief Hot/Cold Gel Sleeve

utilizes temperature holding gel to heat and cool body parts helping to reduce pain from exercise. With its thicker gel and optimal compression, the Mueller EZRelief Hot/Cold Compression Sleeve holds temperature longer, allowing for usage throughout the day, before or after workouts, and during rehabilitation. Intended for reducing swelling and pain, the sleeve can be used on knees, elbows, biceps, wrists, hands, ankles, feet and more.

Ideal for baseball, softball, basketball, football, volleyball, soccer, and more. Pairs well with Mueller Recoverycare Revive Compression Therapy.

Sports Medicine
6618 OSFM
Sportcare
6614 S / 6615 OSFM / 6616 L

Sizing

Measure the largest circumference of the limb for intended application

SMALL 9.0 – 19.3 in / 23.0 – 49.0 cm

OSFM 10.6 – 21.8 in / 26.9 – 55.5 cm

LARGE 11.6 – 23.4 in / 29.4 – 59.5 cm